

PLGF; sFlt; sFlt:PLGF ratio

Pseudonyms: Placental growth factor; soluble FMS like tyrosine kinase-1

The sFlt:PLGF ratio is measured in pregnancy only. These two circulating placentally- derived biomarkers, soluble FMS like Tyrosine kinase 1 (sFlt) and Placental Growth Factor (PLGF), are detectable in the blood of pregnant women and the levels of these markers are altered in pre-eclampsia. The sFlt-1/PLGF ratio, used with standard clinical assessment and subsequent clinical follow-up, is recommended to help rule-out pre-eclampsia in women presenting with suspected pre-eclampsia between 20 weeks and 34 weeks plus 6 days of gestation [NICE DG23].

Pre-eclampsia complicates around 3% of pregnancies and is associated with significant maternal and neonatal morbidity and mortality. Approximately 30% of cases of pre-eclampsia require a preterm delivery to prevent severe maternal complications.

Assessment of suspected pre-eclampsia includes clinical assessment of maternal hypertension, proteinuria, clinical symptoms such as headache, oedema, visual disturbances, foetal growth restriction and measurement of sFlt:PLGF ratio. A high sFlt:PLGF ratio is seen in pre-eclampsia.

General information

Collection container: Serum (with gel separator, 4.9mL brown top Sarstedt tube)

Type and volume of sample: The tubes should be thoroughly mixed before transport to the lab. 1mL whole blood is required as a minimum volume if only sFlt:PLGF ratio is requested.

If referred in from external hospital minimum volume is 0.5ml serum. Samples should be analysed within 48hrs if stored at 2-8°C, stable for up to 6 months stored at -20°C.

Specimen transport/special precautions: Samples should not be taken from patients receiving therapy with high biotin doses (i.e. > 5 mg/day) until at least 8 hours following the last biotin administration.

Apex Test code: PLGFP for ratio (tests not requestable individually)

Laboratory information

Method principle: PLGF and sFlt are analysed on the automated instruments by competitive immunoassay with electrochemiluminescence detection.

Biological reference ranges: PLGF and sFlt both measured in pg/mL, ratio no units

Division of Laboratory Medicine

Biochemistry

Ratio is not reported with reference range as interpretation is based on clinical pathway including clinical maternal and foetal assessment in conjunction with sFlt:PLGF ratio result, not as an isolated finding.

Turnaround times: Results should be available the same working day.

A request can be added on for this test to a sample collected no older than 48hrs.

Clinical information

Factors known to significantly affect the results: None

Clinical decision points: sFlt:PLGF ratio is used in suspected pre-eclampsia not confirmed by clinical assessment <37 weeks. Clinical pathway includes clinical assessment of maternal symptoms, bp, proteinuria, other markers and foetal assessment in conjunction with sFlt:PLGF ratio:

sFlt:PLGF ratio <38	pre-eclampsia excluded for next 7 days
sFlt:PLGF ratio \geq 38 to <85	Intermediate result, increased surveillance
sFlt:PLGF ratio \geq 85	diagnosis of pre-eclampsia/placental disease confirmed

Note a positive ratio is not an isolated indication for delivery <37 weeks.

(Last updated February 2021)