

Division of Laboratory Medicine

Serology

Chlamydia trachomatis antibody

General information

Collection container (including preservatives): 6mL clotted blood tube

Specimen type: venous blood

Collection: 6mL clotted blood tube

Specimen transport: No special needs

Type and volume of sample: 2mL

Special precautions: All samples are suitable for overnight refrigeration only, they must not be stored over a weekend

Laboratory information

Turnaround time: 3 days

Clinical information

Clinical decision points: Not applicable

Factors known to significantly affect the results: Haemolysis

(Last updated December 2014)