

Cytomegalovirus (CMV) IgG avidity

General information

Collection container (including preservatives): 6mL clotted blood tube

Specimen type: Venous Blood

Collection: 6mL blood tube

Specimen transport: No special needs

Minimum volume of sample: 2ml Venous Blood

Special precautions: If processing is delayed, refrigeration is preferable to storage at ambient temperature. Delays of over 48hr are undesirable.

Laboratory information

Turnaround time for provisional result (working days): 5 days

Turnaround time to final result (working days): 7 days

Clinical information

Clinical decision points: Not applicable

Factors known to significantly affect the results: Haemolysis

(Last updated December 2014)