

Epstein Barr Virus (EBNA) - confirmation

General information

Collection container (including preservatives): 6mL clotted blood tube

Specimen type: Venous Blood

Collection: 6mL blood tube

Specimen transport: Ambient temperature via porter, courier, Royal Mail or DX compliant with IATA packing instruction 650

Minimum volume of sample: 2 mL

Special precautions: If processing is delayed, refrigeration is preferable to storage at ambient temperature. Delays of over 48hr are undesirable.

Laboratory information

Measurement units: U/mL

Turnaround time: 3- 4 working days

Clinical information

Clinical decision points: Not applicable

Factors known to significantly affect the results: Haemolysis

(Last updated December 2014)