

## Legionella urinary antigen detection

### General information

**Collection container (including preservatives):** CE marked leak proof container



**Specimen type:** Urine

**Collection:** Urine

**Specimen transport:** No special needs

**Minimum volume of sample:** 1mL urine

**Special precautions:** If processing is delayed, refrigeration is preferable to storage at ambient temperature. Delays of over 48hr are undesirable.

### Laboratory information

**Turnaround time:** Please contact laboratory to arrange urgent test (3 hrs)  
Routine testing is 1 working day

### Clinical information

**Clinical decision points:** Not applicable

**Factors known to significantly affect the results:** Not known

**(Last updated December 2014)**