

Division of Laboratory Medicine

Virology

Lyme IgG, IgM

General information

Collection container (including preservatives):

6mL clotted blood tube
CE marked leak proof container

Specimen type: Venous Blood; CSF can be used for LYME IgG

Collection: 6mL blood tube

Specimen transport: Ambient temperature via porter, courier, Royal Mail or DX compliant with IATA packing instruction 650

Minimum volume of sample: 2 mL

Special precautions: If processing is delayed, refrigeration is preferable to storage at ambient temperature. Delays of over 48hr are undesirable.

Laboratory information

Measurement units: Index

Turnaround time for provisional result (working days): 3 days

Turnaround time to final result (working days): 4 days

Clinical information

Clinical decision points: Not applicable

Factors known to significantly affect the results: Haemolysis

(Last updated December 2014)