

## Division of Laboratory Medicine

Virology

# Pneumococcal urinary antigen detection

## General information

**Collection container (including preservatives):** 10mL urine



**Specimen type:** Urine

**Collection:** 10mL urine

**Specimen transport:** No special needs

**Minimum volume of sample:** 10mL urine

## Laboratory information

**Measurement units:** Not applicable

**Biological reference units:** Not applicable

**Turnaround time for provisional result (working days):** 1 day

**Turnaround time to final result (working days):** 2-3 days

## Clinical information

**Clinical decision points:** Not applicable

**Factors known to significantly affect the results:** None known

**(Last updated December 2014)**