

## Division of Laboratory Medicine

Virology

# Varicella-zoster IgM

## General information

**Collection container (including preservatives):** 6mL clotted blood tube

**Specimen type:** Venous Blood

**Collection:** 6mL blood tube

**Specimen transport:** No special needs

**Minimum volume of sample:** 2 mL

**Special precautions:** None known

## Laboratory information

**Turnaround time for provisional result (working days):** 3 days

**Turnaround time to final result (working days):** 4 days

## Clinical information

**Clinical decision points:** Not applicable

**Factors known to significantly affect the results:** Haemolysis

**(Last updated December 2014)**