

Division of Laboratory Medicine

Immunology

Myelin Associated Glycoprotein (MAG)

Referred to:

Department of Immunology
The Churchill Hospital
Old Road
Headington
Oxford
OX3 7LJ

General information

Specimen transport: At room temperature

Repeat frequency: At significant change of clinical symptoms

Special precautions: None

Laboratory information

Normal reference range: Neg = <1000 BTU (Boulman Time Unit)

Volume and sample type: 4ml serum

Method: ELISA

Turnaround time (calendar days from sample receipt to authorised result): Median - 16, 95th percentile - 34

Clinical information

Indications for the test: IgM monoclonal neuropathy

Factors affecting the test: None

(Last updated October 2023)