

## Division of Laboratory Medicine

Virology

# Hepatitis E IgM

## General information

**Collection container (including preservatives):** 6mL clotted blood tube

**Specimen type:** Venous Blood

**Specimen transport:** Ambient temperature via porter, courier, Royal Mail or DX compliant with IATA packing instruction 650

**Minimum volume of sample:** 2 mL

**Special precautions:** If processing is delayed, refrigeration is preferable to storage at ambient temperature. Delays of over 48hr are undesirable.

## Laboratory information

**Measurement units:** Not applicable

**Turnaround time:** 3-4 working days

## Clinical information

**Clinical decision points:** Not applicable

**Factors known to significantly affect the results:** Highly Lipaemic, icteric and haemolysed samples.

**(Last updated November 2025)**