

## Division of Laboratory Medicine

Immunology

# Amphiphysin

## General Information

**Specimen transport:** At room temperature

**Repeat frequency:** At significant change of clinical symptoms

**Special precautions:** None

## Laboratory information

**Normal reference range:** Negative

**Volume and sample type:** 7ml clotted blood

**Method:** Immuno-Blot

**Turnaround time (calendar days from sample receipt to authorised result):** Median - 16, 95th percentile - 34

## Clinical Information

**Indications for the test:** Paraneoplastic neurological syndromes

**Factors affecting the test:** none

**(Last updated October 2014)**