

## **Division of Laboratory Medicine**

**Immunology** 

# **PLA2R Antibodies**

#### **General information**

This test is referred by us to the Protein Reference Unit at Sheffield. A request for this test on ICE will also generate a label for a TSP7a which is performed in Renal Research Laboratory MFT. Please take two samples and attach a label to each.

For further information see the Protein Reference Unit at Sheffield: https://www.immqas.org.uk/pru.asp?sa=pla2r

Specimen transport: At room temperature

# **Laboratory information**

Normal range: See <a href="https://www.immqas.org.uk/pru.asp?sa=pla2r">https://www.immqas.org.uk/pru.asp?sa=pla2r</a>

**Volume and sample type:** Brown Serum 4.9 mls x2

Method: See <a href="https://www.immgas.org.uk/pru.asp?sa=pla2r">https://www.immgas.org.uk/pru.asp?sa=pla2r</a>

### **Clinical information**

Indications for the test: see <a href="https://www.immgas.org.uk/pru.asp?sa=pla2r">https://www.immgas.org.uk/pru.asp?sa=pla2r</a>