

## Division of Laboratory Medicine

Immunology

### PLA2R Antibodies

**Referred to:**

Department of Immunology  
PO Box 894  
Sheffield  
S5 7YT

#### General information

A request for this test will also generate a label for a TSP7a which is performed in Renal Research Laboratory MFT. Please take two samples and attach a label to each.

For further information see the Protein Reference Unit at Sheffield:

<https://www.immqas.org.uk/pru.asp?sa=pla2r>

**Specimen transport:** At room temperature

#### Laboratory information

**Normal range:** See <https://www.immqas.org.uk/pru.asp?sa=pla2r>

**Volume and sample type:** Brown Serum 4.9 mls x2

**Method:** See <https://www.immqas.org.uk/pru.asp?sa=pla2r>

#### Clinical information

**Indications for the test:** see <https://www.immqas.org.uk/pru.asp?sa=pla2r>

(Last updated January 2025)